

Dia-Stop Cushion

Intended Use

Dia-Stop is a compressive device for puncture sites.

Patient population

It is recommended to use the product on all groups of patients dedicated to haemodialysis treatment with recommended minimum arm circumference equal 15 cm and maximum 33 cm, except for children.

Product description

Dia-Stop is a single use device designed to create haemostasis in the puncture sites established for haemodialysis treatment. Dia-Stop REF 8065-Cn consists of single use sterile compression cushion. Combined with reusable plastic armlet with velcro attachment (REF 8065-A) it creates Dia-Stop compressive device. It has a unique bridge design that ensures that no pressure is applied to a substantial part of the circumference of the arm. This eliminates the risk of stasis/interruption of the blood flow. The pressure cushion of the Dia-Stop is made of soft foam materials and includes a pressure indicator. Compression of the dark cushion to less than half thickness indicates a suitable pressure is applied on an average patient. This minimizes the risk of causing damage to the fistula.

Precautions

- Read the manufacturer's instruction prior to use.
- For single use only. Do not reuse, reprocess or re-sterilize this device. Any reprocessing may impede the functions of the device. Reusing single use devices may also increase the risk of cross contamination. Attempts to clean the device results in risk of device malfunction.
- Carefully inspect the packaging for any damage prior to use. Do not use if the packaging is damaged. Do not use past expiration date. Do not use if the device is exposed to non-sterile surface before procedure.
- After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.
- Dia-Stop should be used by trained and/or experienced person.

Instruction for Use

In order to achieve haemostasis in the puncture site please follow below steps:

1. Take the Dia-Stop armlet out of the pouch (Figure 1).
2. Wipe the armlet thoroughly using 70 % v/v Isopropyl Alcohol wipes (Figure 2).
3. Take the Dia-Stop cushion out of the pouch (Figure 3).
4. Attach the cushion by the Velcro hook onto the armlet (Figure 4).
5. Place the armlet with cushion around arm and pull the notches (Figure 5).
6. Take the needle out of the puncture site (1) and applied pressure (2) onto the compression cushion (Figure 6).

7. Adjust the notches in order to obtain compression of the dark cushion to less than half its initial thickness (Figure 7).
8. After adjustment of the notches to correct size, hook the notch in the small hole (Figure 8).
9. Dia-Stop application is finished.

In order to remove and then reuse the armlet:

1. Pull up the rectangle-shaped end of the armlet (Figure a).
2. Remove the used compression cushion and dispose (Figure b).
3. Wipe the armlet thoroughly using 70 % v/v Isopropyl Alcohol wipes (Figure c).
4. Take Dia-Stop cushion out of the pouch (Figure d).
5. Follow again steps from 4 to 9 from above instruction.

It is recommended to use one armlet up to 10 times.

Reporting

Any serious incident that has occurred in relation to the device should be reported to the LiNA Medical ApS and the competent regulatory authority of the country in which the user and/or patient is established.

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