



Innovation in Gynaecology



LiNA EasyBag™

The cost effective specimen retrieval system





LiNA EasyBag™



safe

Strong/tear resistant pouch

The pouch of the LiNA EasyBag is made of strong polyurethane plastic. The LiNA EasyBag is available in two different sizes; 75 x 150mm (250ml) and 125 x 200mm (750ml).

functional

Easy to use

The LiNA EasyBag is enclosed inside a 10mm holster for easy insertion into a standard trocar. After insertion unfold the LiNA EasyBag with a grasper and fill in the specimen. Remove the LiNA EasyBag through the incision, after removal of the trocar.

versatile

Use the LiNA EasyBag for the following procedures:

- Appendices
 - Fallopian Tubes
 - Gall bladder and stones
 - Section of bowel
 - Omentum
 - Spleen
 - Blood Clots
 - Lymph nodes
 - Lung wedges
 - Fibroid Myomas
 - Ovaries
 - Kidneys
- and other similar structures

Ordering Information:

LiNA EasyBag **CE** 0459

REF.	Product	Packing
EB-75	LiNA EasyBag - small, 75 x 150mm, volume 250ml. Single use.	Box of 20, sterile
EB-125	LiNA EasyBag - large, 125 x 200mm, volume 750ml. Single use.	Box of 20, sterile

