



Innovation in Gynaecology



LiNA EndoHook™

Single use endoscopic monopolar hook electrodes for laparoscopic procedures



- *Available with 3 different tips*
- *Choose between finger- or footswitch instruments*

safe

- Ceramic tip insulation designed to lead the energy exactly where it is needed

versatile

Available in different versions


- Fingerswitch or footswitch
- Length 320 mm or 400 mm
- Needle tip, L-shaped tip or Spatula tip electrode
- 4 mm male connector or 3 pin connector
- Adaptors for different generators



functional

- Lightweight and slim design
- Different tip designs for precise dissection, cutting and coagulation in a variety of surgical applications

LiNA Tip cleaner
For cleaning of electrodes.



LiNA EndoHook™ **CE 0459**

Ordering information:

REF.	Product	Packing
308-L (L-shaped tip) 308-N (Needle tip) 308-S (Spatula tip)	LiNA EndoHook - 320 mm instrument – Foot controlled instrument with 4 mm male connector in the handle, excl. cable – Single use. Available with 3 different tips.	Box of 10, sterile.
308-L-400 (L-shaped tip)	LiNA EndoHook - As above but with extra long 400 mm instrument.	Box of 10, sterile.
572 (for 308)	LiNA HF generator Adaptor – reusable – 8 mm connector for foot controlled endo hooks. (for Valleylab generators).	Box of 1.
REC-617 (for 308)	Monopolar Silicone Cable - w. 4 mm connector, 3 m, reusable.	Box of 1.
908-L (L-shaped tip) 908-N (Needle tip) 908-S (Spatula tip)	LiNA EndoHook – 320 mm instrument – Finger controlled instrument with 3 m integrated cable and 3 pin standard generator connector – Single use. Available with 3 different tips.	Box of 10, sterile.
571-ACC (for 908) 571-AE 571-B 571-E 571-L 571-M 571-N 571-S1 571-S2	LiNA HF generator Adaptor – reusable – For Erbe ACC generator LiNA HF generator Adaptor – reusable – For Aesculab generator LiNA HF generator Adaptor – reusable – For Bovie generator LiNA HF generator Adaptor – reusable – For Old Erbe generator (Erbotom T) LiNA HF generator Adaptor – reusable – For Lamidey generator LiNA HF generator Adaptor – reusable – For Martin generator LiNA HF generator Adaptor – reusable – For Neomed generator LiNA HF generator Adaptor – reusable – For Siemens type 1 generator LiNA HF generator Adaptor – reusable – For Siemens type 2 generator	Box of 1.
580	LiNA Tip Cleaner - Single use - Self Adhesive. 45x50 mm with X-ray detectable thread.	Box of 25, sterile