

# M/L-10

MULTI-FIRE LAPAROSCOPIC CLIP APPLIER

Get a handle on  
cost-conscious  
reposability.



**MICROLINE**  
SURGICAL

*The Reposables Company*

# M/L-10

MULTI-FIRE LAPAROSCOPIC CLIP APPLIER


The M/L-10 Clip Applier is a reusable multi-fire device with a high-quality reusable handle. It employs disposable cartridges that contain 19 implant-grade titanium clips, the last two of which are color-coded.



Procedure Kit

## Ligating Clip Cartridge

- Unique cross pattern design
- 2 to 5 mm vessel/cystic duct ligation capability
- Surgical-grade titanium
- Tips-first closure sequence for maximum ligating capability and dependability
- Consistent, uniform closure



### ORDERING INFORMATION:

M/L-10 Multi-Fire Clip Applier	REF 1002
M/L Clip Cartridge (19 clips)	REF 1112
Spare Outer Tube	REF 1201

### Cost Saving Procedure Kit:

Clip Cartridge with 3142 Endo Cut Scissor Tip	REF 1242
---	----------

**MICROLINE**  
SURGICAL

50 Dunham Road | Beverly, MA 01915 USA  
Phone 978-922-9810 | Fax 978-922-9209 | [www.microlinesurgical.com](http://www.microlinesurgical.com)

© 2015 Microline Surgical. Microline is a trademark and/or registered trademark of Microline Surgical.  
P/N: 09-39-00054-00-US Rev E, 6/15

**Important information:** Prior to use, refer to Instructions for Use supplied with this device for indications, contraindications, adverse events, suggested procedures, warnings and precautions for full product information.

**Caution:** Federal (USA) law restricts this device to sale by or on the order of a Physician.