

Large Sampling Capability
Leads to Excellent Detection Ability

ANSO[®] TRIIPLE NEEDLE BRUSH

SINGLE-USE
ⓧ





We listen, Innovate, and Create to save lives

“Three brushes, triple
specimen collection”

ANSO® TRIPLE NEEDLE BRUSH

CANCER DIAGNOSTIC



| FEATURES | BENEFITS |
|--|---|
| Unique design | Three needle brushes vs one needle brush |
| Larger Sampling of cells collection in a single pass | Three brushes with needle, vs one needle brush, triples the specimen collection |
| Better Accuracy | Larger sampled area, produces larger cell collection. |
| Excellent Suction Capability | Triple area for suctioning |

PRODUCT SPECIFICATIONS ANSO® TRIPLE NEEDLE BRUSH

| PRODUCT CODE | CATETER DIAMETER MM | CATETER LENGTH CM | BRUSH DIAMETER LENGTH MM | PACKAGING |
|--------------|---------------------|-------------------|--------------------------|-----------|
| HIC-CB-27-10 | 1.8 | 120 | 0.87 / 15 | BOX 10 |
| HIC-CB-27-5 | 1.8 | 120 | 0.87 / 15 | BOX 5 |

ANSO® is a registered Trademark of Horizons International Corp.

Distributed by:



We listen, Innovate, and Create to save lives

Horizons International Corp.

P.O. Box 213-3006 Zona Franca Metro,
Barreal de Heredia, Costa Rica

Phone: (506) 2261-9114 / Fax: (506) 2261-9115

FOR ORDERS IN THE U.S. PLEASE CALL TOLL FREE 866-245-7697

website: www.horizonscorp.com / **email:** info@horizonscorp.com