


# COBRA™

The Power of Two in One  
Rotatable Snare / Retrieval Device



360° Rotation





Horizons International Corp.

*We listen, Innovate, and Create to save lives*

# COBRA™

## The Power of Two in One

### Rotatable Snare / Retrieval Device



| FEATURES                                    | BENEFITS   |
|---|--|
| The Power of Two in One                     | Unique Patented design, half snare and half Rerievel Net, to cut and collect polyps in one single pass |
| 360° Rotation (1:1 ratio in all directions) | Designed to facilitate cutting , capturing and removing of polyps, in one single pass                  |
| Deep, Three - Dimensional Net               | Engineered to capture and remove multiple polyps, in a single pass                                     |
| Durability: Interwoven Net Design           | Net is attached to the outer wire using a interwoven design to help maintain the integrity of the net  |
| Short throw handle                          | Ergonomically design for user comfort  |

#### PRODUCT SPECIFICATIONS

### COBRA™ Rotatable Snare / Retrieval Device

| PRODUCT CODE   | CATETER DIAMETER MM | SIZE MM | LENGTH CM | TYPE      | PACKAGING |
|----------------|---------------------|---------|-----------|-----------|-----------|
| HIC-SODT-13-5  | 2.5                 | 13      | 230       | OVAL      | BOX 5     |
| HIC-SODT-13-10 | 2.5                 | 13      | 230       | OVAL      | BOX 10    |
| HIC-SODT-24-5  | 2.5                 | 25      | 230       | OVAL      | BOX 5     |
| HIC-SODT-24-10 | 2.5                 | 25      | 230       | OVAL      | BOX 10    |
| HIC-SODT-25-5  | 2.5                 | 32      | 230       | OVAL      | BOX 5     |
| HIC-SODT-25-10 | 2.5                 | 32      | 230       | OVAL      | BOX 10    |
| HIC-SHDT-13-5  | 2.5                 | 13      | 230       | HEXAGONAL | BOX 5     |
| HIC-SHDT-13-10 | 2.5                 | 13      | 230       | HEXAGONAL | BOX 10    |
| HIC-SHDT-24-5  | 2.5                 | 25      | 230       | HEXAGONAL | BOX 5     |
| HIC-SHDT-24-10 | 2.5                 | 25      | 230       | HEXAGONAL | BOX 10    |
| HIC-SHDT-25- 5 | 2.5                 | 32      | 230       | HEXAGONAL | BOX 5     |
| HIC-SHDT-25-10 | 2.5                 | 32      | 230       | HEXAGONAL | BOX 10    |

COBRA™ SNARE is a Trademark of Horizons International Corp

Protected by one or more of the following US Patents: 9,101,342; 8,858,567. One patent pending.

Caution: Federal (USA) law restricts this device to be sold by, or on the order of a physician.

Distributed by:



*We listen, Innovate, and Create to save lives*

#### Horizons International Corp.

P.O. Box 213-3006 Zona Franca Metro,

Barreal de Heredia, Costa Rica

Phone: (506) 2261-9114 / Fax: (506) 2261-9115

**FOR ORDERS IN THE U.S. PLEASE CALL TOLL FREE 866-245-7697**

**website:** [www.horizonscorp.com](http://www.horizonscorp.com) / **email:** [info@horizonscorp.com](mailto:info@horizonscorp.com)